



Sarcopenia Assessment

using DXA Technology from GE Healthcare

Healthy Aging – It's Vital



What is sarcopenia?

Sarcopenia – muscle loss with aging

Sarcopenia is a disease associated with the aging process. Loss of muscle mass and strength affects balance, gait, and the ability to perform daily living tasks. Sarcopenia most commonly is seen in inactive people but can also affect individuals who remain physically active throughout their lives.

Although a sedentary lifestyle contributes to this disease, it's not the only factor. As we age, hormone levels change, protein requirements alter, motor neurons die, and we tend to become more sedentary.¹

In recent studies, it has been found that worsening sarcopenia followed trends in losses of muscle strength as well as impairment of daily functioning.² In one such study, it was found that the prevalence of sarcopenia increased dramatically with age from 4% of men and 3% of women aged 70-75 to 16% of men and 13% of women aged 85 or older.³

When sarcopenia is coupled with other diseases associated with aging, its effects can be even more pronounced. Loss of muscle mass and strength is a significant risk factor for disability in the aging population.⁴ When patients suffer from both sarcopenia and osteoporosis, the risk of falling and subsequent fracture incidence is higher.⁵ Therefore, treating sarcopenia will in turn help to lessen its burden on co-existing diseases.




Illustration of sarcopenia induced muscle loss



DXA bone densitometer systems with enCORE v17 sarcopenia feature

The optional sarcopenia software in enCORE version 17 calculates values based on published definitions and thresholds using measured **appendicular lean mass** in combination with patient demographics and entered values of **muscle strength** and **physical performance**. These values may be useful to health care professionals in their management of sarcopenia.


Sarcopenia Calculator uses Appendicular Lean Mass (ALM) from a DXA Scan



Choose from different sarcopenia definitions

The sarcopenia feature allows you to choose from one of the five standards:


- **BAUM** – Baumgartner⁶
- **AWGS** – Asian Working Group for Sarcopenia⁷
- **EWGSOP** – European Working Group on Sarcopenia in Older People⁸
- **FNIH** - Foundation for National Institutes of Health Sarcopenia Project⁹
- **IWGS** – International Working Group on Sarcopenia¹⁰



Possible Definitions of Sarcopenia in enCORE 17

View an existing sarcopenia definition

You could also view any of the existing sarcopenia definitions.



Sarcopenia Definitions

Enable

| | |
|---------------|---|
| AWGS | Asian Working Group for Sarcopenia |
| BAUM | Baumgartner |
| EWGSOP | European Working Group on Sarcopenia in Older People |
| FNIH | Foundation for the National Institutes of Health Sarcopenia Project |
| IWGS | International Working Group on Sarcopenia |

Advanced

OK
Cancel
New...
View...

View existing Sarcopenia definition

Sarcopenia Definition

Abbreviation and Name: EWGSOP European Working Group on Sarcopenia in Older People

Female:

Minimal Age: 65

| Sev. | Classification | Formula |
|------|-------------------|--|
| 1 | Presarcopenia | $(ALM/Ht^2 \leq 5.45 \text{ kg/m}^2)$ |
| 2 | Sarcopenia | $((\text{Gait Speed} < 0.80 \text{ m/s}) \text{ Or } (\text{Grip Strength} < 20.0 \text{ kg})) \text{ And } (ALM/Ht^2 \leq 5.45 \text{ kg/m}^2)$ |
| 3 | Severe Sarcopenia | $(\text{Gait Speed} < 0.80 \text{ m/s}) \text{ And } (\text{Grip Strength} < 20.0 \text{ kg}) \text{ And } (ALM/Ht^2 \leq 5.45 \text{ kg/m}^2)$ |

Male:

Minimal Age: 65

| Sev. | Classification | Formula |
|------|-------------------|--|
| 1 | Presarcopenia | $(ALM/Ht^2 \leq 7.26 \text{ kg/m}^2)$ |
| 2 | Sarcopenia | $((\text{Gait Speed} < 0.80 \text{ m/s}) \text{ Or } (\text{Grip Strength} < 30.0 \text{ kg})) \text{ And } (ALM/Ht^2 \leq 7.26 \text{ kg/m}^2)$ |
| 3 | Severe Sarcopenia | $(\text{Gait Speed} < 0.80 \text{ m/s}) \text{ And } (\text{Grip Strength} < 30.0 \text{ kg}) \text{ And } (ALM/Ht^2 \leq 7.26 \text{ kg/m}^2)$ |

OK Cancel



Create new sarcopenia definitions

You also have the ability to create new sarcopenia definitions by using patient measurement parameters such as ALM, BMI, height, gait speed and grip strength.

- **Classification field:** Enter a classification name
- **Logic drop-down menus:** Provide the classification measurement device “and/or” logic. Select Or() and And() to provide parentheses to indicate precedence
- **Measurement:** Select measurement devices used in determining sarcopenia classifications
- **Operator selections:** Provide a math operation for comparison to value
- **Value fields:** Enter a comparison measurement value

Sarcopenia Definition

Abbreviation and Name: USER | New Sarcopenia Definition

Female: Minimal Age | 65

Male: Minimal Age | 65

Classification New

Classification: New Sarcopenia Definition

Logic: And

Measurement: Gait Speed

Operator: < (is less than)

Value: 0.01 m/s

NOTE: Classifications must be ordered from least to highest severity (Sev.)

Classification New

Classification: New Sarcopenia Definition


| Logic: | Measurement: | Operator: | Value: |
|--------|---------------------|-------------------------------|------------------------|
| And | Gait Speed | < (is less than) | 0.80 m/s |
| And | Grip Strength | < (is less than) | 20.0 kg |
| Or | ALM/HT ² | <= (is less than or equal to) | 5.54 kg/m ² |

OK Cancel



Calculate sarcopenia results

After you have selected a default sarcopenia definition, you can use the sarcopenia feature to calculate whether a patient is sarcopenic after entering values of muscle strength and physical performance.



Sarcopenia Calculator



References:

1. International Osteoporosis Foundation - <https://www.iofbonehealth.org/what-sarcopenia>
2. Mithal A, Bonjour J-P, Boonen S, Burckhardt P, Degens H, El Hajj Fuleihan G, Josse R, Lips P, Morales Torres J, Rizzoli R, Yoshimura N, Wahl D.A., Cooper C, Dawson-Hughes B(2011) Impact of nutrition on muscle strength and performance in older adults. Osteoporosis International (in press)
3. Castillo EM, Goodman-Gruen D, Kritz-Silverstein D, Morton DJ, Wingard DL, Barrett-Connor E (2003) Sarcopenia in elderly men and women - The Rancho Bernardo Study. Am J Prev Med 25: 226-231
4. Volpi, E, Nazemi R, Fujita S, (2004) Muscle tissue changes with aging. Curr Opin Nutr Metab Care July, 7(4):405-410
5. Sarcopenia: European consensus on definition and diagnosis, Report of the European Working Group on Sarcopenia in Older People," Age and Ageing Advance Access originally published online on April 13, 2010, Age and Ageing 2010 39(4):412-423
6. Baumgartner RN, Koehler KM, Gallagher D et al (1998) Epidemiology of sarcopenia among the elderly in New Mexico. Am J Epidemiol 147(8):755-763
7. Chen LK, Liu LK, Woo J, Assantachai P, Auyeung TW, Bahyah KS, et al. Sarcopenia in Asia: consensus report of the Asian Working Group for Sarcopenia. JAMDA. 2014;15:95-101
8. Cruz-Jentoft AJ, Baeyens JP, Bauer JM et al (2010) Sarcopenia: European consensus on definition and diagnosis: report of the European Working Group on Sarcopenia in Older People. Age and Ageing 39(4):412-423
9. Studenski SA, Peters KW, Alley DE et al (2014) The FNIH sarcopenia project: rationale, study description, conference recommendations, and final estimates. J Gerontol A Biol Sci Med Sci 69(5):547-558
10. Fielding RA, Vellas B, EvansWJ et al (2011) Sarcopenia: an undiagnosed condition in older adults. Current consensus definition: prevalence, etiology, and consequences. International working group on sarcopenia. J Am Med Dir Assoc 12(4):249-256



© 2017 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, and the GE Monogram are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

April 2017
JB48544XX