



Bone & Metabolic Health

# Advanced Analytics

Advanced insights into bone  
& body composition

[gehealthcare.com](http://gehealthcare.com)

# enCORE v18 Advanced Analytics

The latest upgrade to our enCORE software elevates our world-class DXA family of products to a whole new level

Advanced Analytics allows you to create custom equations using existing data from your DXA system, set user-defined classification thresholds based on these equations, and establish actionable patient goals based on evaluation of the results. Later, these custom ratios can be printed out in a custom-built report to share the results with end users.

# Hundreds of Bone and Body Composition Parameters are available

Using Advanced Analytics, Clinicians, Sports Meds & Body Composition researchers can analyze DXA bone & body composition data using sophisticated techniques to discover deeper insights, create custom ratios, make predictions, and have metrics pinned on an in-built dashboard in enCORE software to track changes over time.

Metrics and trending can be performed on retrospective patient or athlete data too.

Advanced Analytics empowers researchers and clinicians to create a “metric” based on more than 200 DXA bone and body composition parameters captured during a DXA scan allowing creation of custom ratios and simple or complex mathematical equations.

Custom metrics can be used for research or in clinical practice to gain deeper insights into the changes in bone mineral density and body composition of various regions of interest over time.

## MORE THAN 200 BONE AND BODY COMPOSITION PARAMETERS FOR ANALYSIS.

Custom Regions of Interest (ROI) can also be used with Advanced Analytics feature.

ADVANCED ANALYTICS  
TOOLS POWER DEEPER  
INSIGHTS INTO BONE AND  
BODY COMPOSITION

Sample Metric = VAT Mass/Total Fat Mass

BODY COMPOSITION ANALYTICS			BONE ANALYTICS – VARIOUS SKELETAL SITES		
DXA Parameters BMC Fat Free Mass Fat Mass Lean Mass Region %Fat Tissue %Fat Tissue Mass Total Mass	Regions Left Arm Right Arm Android Gynoid VAT SAT Total Body TBLH much more	Composition		AP Spine	Regions L1 L2 L3 L4 L1-L2 L1-L4 etc.
	Regions Left Arm Right Arm Android Gynoid VAT SAT Total Body TBLH much more				DXA Parameters Area BMD BMC BMD %AM BMD %YA BMD T-Score BMD Z-Score
DXA Parameters Area BMC BMD BMD %AM BMD %YA BMD T-Score BMD Z-Score	Regions Left Arm Right Arm Android Gynoid VAT SAT Total Body TBLH much more	Densitometry		Femur	Regions Neck Trochanter Shaft Total etc.
	Regions Left Arm Right Arm Android Gynoid VAT SAT Total Body TBLH much more				DXA Parameters AFF Beaking Index Area BMD BMC BMD %AM BMD %YA BMD T-Score BMD Z-Score

Note: Some parameters may require purchase of additional features.



# Advanced Insights

Clinicians, Sports Medicine & Body Composition researchers can analyze DXA bone & body composition data using sophisticated techniques to discover deeper insights, create custom ratios, make predictions, have results pinned on a dashboard, track changes of custom ratios over time, or generate clinical recommendations. Customer ratios can be run retrospectively too.

## ADVANCED INSIGHTS INTO BONE HEALTH

- Radiologists
- Rheumatologist
- Orthopedic
- Oncologist
- OB/GYN



- Endocrinologists
- Sports Med
- Nutritionists & Dietitians

## ADVANCED INSIGHTS INTO BODY COMPOSITION

# Create new Metrics

Create simple or complex mathematical formulas using parameters from your bone and body composition exams. Mathematical operators and parentheses have been provided to create complex mathematical equations.

Equation:  
**Fat Mass Index (FMI) (kg/m<sup>2</sup>) = Total Fat Mass (Kg) / Patient Height [m<sup>2</sup>]**

Math Operators

+

-

\*

/

Add

Remove

(

)

(

)

Move

Up

Down

Easily create your own mathematical formulas and test your equation before it is finalized.

# Define Thresholds to drive results from your metrics

Classification Thresholds

Female:

Classification

Underweight

18.5

Normal

25.0

Overweight

30.0

Obesity

Male:

Classification

Underweight

18.5

Normal

25.0


Overweight

30.0

Obesity

Define custom thresholds to enable automatic classifications and evaluation parameters.

# Advanced Analytics Dashboard




Each created metric will be visible in the Analytics Dashboard. Newly created analytics can also be applied retrospectively to existing files.




## Trend Metrics over time

Changes in metrics can be tracked and monitored over time when the patient or athlete comes back for a DXA scan.



## Your Metrics in Custom-Built Reports

Visual results and charts in custom reports that are easy to share. You can set thresholds and classifications to better communicate to patients and clients about changes over time in their bone and body composition.





© 2019 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and DXAVision are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

September 2019  
JB70812XX