

High Purity Aluminum HVL

Enhance the accuracy of your measures.

- Ultra-thin, high purity aluminum
- Ideal for mammographic, fluoroscopic and other radiographic systems



Gammex's ultra-thin and high purity aluminum sheets provide more accurate measures of half value layer aluminum (HVL Al) than other less pure aluminum alloys.

The High Purity Aluminum HVL Set provides the needed purity for accurate MQSA and other HVL compliance testing. The set consists of six 10.0 x 10.0 x 0.1 cm sheets of 99.99% pure aluminum. The set also includes a plastic storage case to help protect and maintain flatness of the thin filters.

Specifications

Construction:	99.99% Pure Aluminum
Quantity:	6
Length:	10.0 cm (4.0 in)
Width:	10.0 cm (4.0 in)
Thickness:	0.1 mm
Weight:	0.09 kg (0.2 lbs)