


Scatter Radiation Diagrams


Isodose Diagram: Prodigy, Prodigy Advance, Prodigy Primo, Prodigy Forma, Prodigy Pro Full-Size Table

The following display isodose diagrams of the Prodigy series full size scanner scatter radiation. The measurements were taken with a Victoreen 470A. The beam was attenuated through a 20 cm water phantom.

Isodose diagram: Foot-end view, Prodigy full-size table



Isodose diagram: Front view, Prodigy full-size table



Legend

— Scanner was run at 3.00 mA.