



CUSTOMER PROVIDED:

- RECEPTACLES:
- OUTLET "A": Outlet, Shall be for the Bone Density Machine, dedicated circuit, 120 volts 20 amps. NEMA 5-20R Dedicated (AS PER LOCAL CODE)
 - OUTLET "B": Network drop and telecommunication outlet and Power Outlet for Computer & Printer. 120 volts 20 amps.




 DOUBLE RECEPTACLE OUTLET

 COMPUTER NETWORK PORT
(100Mbps / 1 Gbps) AND
PHONE RECEPTACLE

iDXA – ROOM LAYOUT
Optimum Sample 14' X 14' Room

Aymes Medical
Date Prepared: Sample Drawing



BMD ROOM

ALL DIMENSIONS = IMPERIAL / [METRIC]